
Bias-UQ-PCR

Release v0.0.1

Jan 12, 2023

CONTENTS:

1	Get Data	1
2	Wells	3
3	Amplification	5
4	Molar Fluorescence	9
5	Kinetic PCR	11
6	Plotting Figures	13
7	Indices and tables	15
	Python Module Index	17

GET DATA

`src.get_data.cycle_to_index(cycle: int) → int`

Parameters `cycle` (*int*) – cycle number *i*, 0 to *n*

Returns index in matrix or array

Return type `int`

`src.get_data.file_to_numpy(file_abs_path: str) → numpy.ndarray`

Parameters `file_abs_path` (*str*) – path to file name

Notes

All fluorescence values are divided by 10^6 before any analysis

Returns *n* by *m* matrix of fluorescence, **F**

Return type `np.ndarray`

`src.get_data.index_to_cycle(index: int) → int`

Parameters `index` (*int*) – index of cycle (0 to *n*-1)

Returns cycle number (1 to *n*)

Return type `int`

WELLS

`src.wells.column_row_to_well (ix: int, iy: int) → str`
Convert indices to well

Parameters

- **ix** (*int*) – x (column) index
- **iy** (*int*) – y (row) index

Returns name of well

Return type str

`src.wells.number_to_column (number: int) → int`
Get column index associated with well number

Parameters **number** (*int*) – well number (0 to 95)

Returns column (0 to 11)

Return type int

`src.wells.number_to_row (number: int) → int`
Get the row number associated with a well number

Parameters **number** (*int*) – well number

Returns row number (0 to 7)

Return type int

`src.wells.number_to_well (number: int) → str`
Get well name from number

Parameters **number** (*int*) – well number (0 to 95)

Returns well name (A1 to H12)

Return type str

`src.wells.well_to_column (well: str) → int`
Convert well name to column

Parameters **well** (*str*) – name of well

Returns **ix** – x-index for well (row)

Return type int

`src.wells.well_to_column_row (well: str) → Tuple[int]`
Convert well name to column, row

Parameters `well` (*str*) – name of well

Returns column index, row index

Return type tuple(int, int)

`src.wells.well_to_number` (*well: str*) → int

Returns number of well

Parameters `well` (*str*) – name of well (e.g., "A1")

Returns well number (0 to 95)

Return type int

`src.wells.well_to_row` (*well: str*) → int

Well to y index

Parameters `well` (*str*) – well name

Returns `iy` – y-index of well (row)

Return type int

AMPLIFICATION

```
class src.amplification.Amplification (R: float, pbar: float, E_U0: numpy.ndarray, V_U0:  
numpy.ndarray, farray=<built-in function array>, sqrt2=1.4142135623730951)
```

Parameters

- **l1** (*float*) – Largest eigenvalue of **A**, λ_1
- **l2** (*float*) – Smallest eigenvalue of **A**, λ_2
- **x1** (*np.ndarray*) – Right eigenvector of **A** corresponding to λ_1 , \mathbf{x}_1
- **x2** (*np.ndarray*) – Right eigenvector of **A** corresponding to λ_2 , \mathbf{x}_2
- **z1** (*np.ndarray*) – Left eigenvector of **A** corresponding to λ_1 , \mathbf{z}_1
- **z2** (*np.ndarray*) – Left eigenvector of **A** corresponding to λ_2 , \mathbf{z}_2
- **K1** (*np.ndarray*) – Matrix **K**₁ defined in text to calculate variance
- **K2** (*np.ndarray*) – Matrix **K**₂ defined in text to calculate variance

```
__init__ (R: float, pbar: float, E_U0: numpy.ndarray, V_U0: numpy.ndarray, farray=<built-in func-  
tion array>, sqrt2=1.4142135623730951) → None
```

Parameters

- **R** (*float*) – ratio of amplification probabilities, R
- **pbar** (*float*) – geometric mean of amplification probabilities, \bar{p}
- **E_U0** (*np.ndarray*) – initial expected values of both strands $\mathbb{E}[\mathbf{U}_0]$, 2 by 1 matrix
- **V_U0** (*np.ndarray*) – initial variances of both strands $\text{Var}[\mathbf{U}_0]$, 2 by 2 matrix
- **farray** (*callable, optional*) – function to convert to array type, defaults to np.array
- **sqrt2** – square root of 2, defaults to float calculated by np.sqrt(2)

```
get_A() → numpy.ndarray
```

Returns matrix **A** calculated by decomposition (10)

Return type np.ndarray

```
get_Atoi (i: int) → numpy.ndarray
```

Parameters **i** (*int*) – cycle number

Returns matrix **A**^{*i*} calculated by decomposition (10)

Return type np.ndarray

get_Exi_over_EYi (*cycles: numpy.array*) \rightarrow *numpy.array*

Calculate $\mathbb{E}[X_i] / \mathbb{E}[Y_i]$ for a variety of cycles i

Parameters **cycles** (*np.array*) – cycles (integers) to calculate for each i

Returns $\mathbb{E}[X_i] / \mathbb{E}[Y_i]$

Return type *np.array*

get_E_Ui (*i: int*) \rightarrow *numpy.ndarray*

Parameters **i** (*int*) – cycle number

Returns matrix $\mathbb{E}[\mathbf{U}_i] = \mathbf{A}^i \mathbb{E}[\mathbf{U}_0]$

Return type *np.ndarray*

get_V_Ui (*i: int*) \rightarrow *numpy.ndarray*

Parameters **i** (*int*) – cycle number

Returns complicated expression to calculate $\text{Var}[\mathbf{U}_i]$

Return type *np.ndarray*

get_nu ()

Returns Computes ν , see Equation (22).

Return type *float*

get_x1z1TK1z1x1T () \rightarrow *numpy.ndarray*

Returns computes $\mathbf{x}_1 \mathbf{z}_1^\top \mathbf{K}_1 \mathbf{z}_1 \mathbf{x}_1^\top$

Return type *np.ndarray*

get_x1z1TK2z1x1T ()

Returns computes $\mathbf{x}_1 \mathbf{z}_1^\top \mathbf{K}_2 \mathbf{z}_1 \mathbf{x}_1^\top$

Return type *np.ndarray*

`src.amplification.eval_E_D0` (*E_U0: numpy.ndarray, R: float*)

Calculates $\mathbb{E}[D_0]$ from $\mathbb{E}[\mathbf{U}_0]$ and R

Parameters

- **E_U0** (*np.ndarray*) – Initial expected value vector
- **R** (*float*) – square root of ratio of amplification probabilities

Returns $\mathbb{E}[D_0] := \mathbb{E}[X_i - RY_i]$

Return type *float*

`src.amplification.eval_E_S0` (*E_U0: numpy.ndarray, R: float*)

Calculates $\mathbb{E}[S_0]$ from $\mathbb{E}[\mathbf{U}_0]$ and R

Parameters

- **E_U0** (*np.ndarray*) – Initial expected value vector
- **R** (*float*) – square root of ratio of amplification probabilities

Returns $\mathbb{E}[S_0] := \mathbb{E}[X_i + RY_i]$

Return type *float*

`src.amplification.eval_R(p_fr: float, p_rf: float)`

Calculate R from p_{rf} and p_{fr}

Parameters

- **p_fr** (*float*) – probability of forward to reverse amplification, p_{fr}
- **p_rf** (*float*) – probability of reverse to forward amplification, p_{rf}

Returns $R = \sqrt{\frac{p_{rf}}{p_{fr}}}$

Return type *float*

`src.amplification.get_pfr_prf(R: float, pbar: float)`

Get p_{fr} and p_{rf} from R and \bar{p}

Parameters

- **R** (*float*) – square root of ratio of probabilities, R
- **pbar** (*float*) – geometric mean of probabilities, \bar{p}

Returns

Return type *tuple(p_fr, p_rf)*

MOLAR FLUORESCENCE

```
class src.molar_fluorescence.MolarFluorescence (C: numpy.array, files: List[str], name:  
                                              str)
```

Parameters

- **f** (*np.ndarray*) – molar fluorescences for each cycle/well **f**, **n** by number of wells matrix
Determined in units of fluorescence divided by (pmol/L)
- **df** (*np.ndarray*) – standard deviation in molar fluorescences σ , **n** by number of wells
matrix Determined in units of fluorescence divided by (pmol/L)
- **q** (*int*) – number of plates (different concentrations)
- **n** (*int*) – number of cycles
- **m** (*int*) – number of wells
- **F** (*np.ndarray*) – raw fluorescence data (scaled by 10^6 as described in `get_data.py`).
Tensor of dimensions $n \times m \times q$

```
__init__ (C: numpy.array, files: List[str], name: str)
```

Parameters

- **C** (*np.array*) – concentrations of reporter in pmol/L
- **files** (*typing.List[str]*) – list of file names (absolute paths)
- **name** (*str*) – name of reporter (e.g., FAM, Probe)

```
calculate ()
```

Perform calculations of **f** and standard deviation σ

KINETIC PCR

```
class src.kinetic_PCR.HydrolysisProbes(C: float, Vol: float, f_plus: numpy.ndarray,
                                       f_minus: numpy.ndarray, R: float, pbar: float, E_U0:
                                       numpy.ndarray, V_U0: numpy.ndarray, **kwargs)
```

Specific subclass for hydrolysis probes

Parameters

- **n** (*int*) – number of cycles
- **m** (*int*) – number of wells
- **d** (*np.ndarray*) – array of incremental increases in fluorescences, n by m
- **b** (*np.ndarray*) – array of background signals, n by m
- **E_F** (*np.ndarray*) – expected value of fluorescence, n by m, $\mathbb{E}[\mathbf{F}]$
- **V_F** (*np.ndarray*) – variance of fluorescence, n by m, $\text{Var}[\mathbf{F}]$
- **E_DX** (*np.ndarray*) – expected value of change in X, length n, $\mathbb{E}[\Delta X_i]$
- **V_DX** (*np.ndarray*) – variance value of change in X, length n, $\text{Var}[\Delta X_i]$
- **cycles** (*np.ndarray*) – cycle numbers, 1 to n

```
__init__(C: float, Vol: float, f_plus: numpy.ndarray, f_minus: numpy.ndarray, R: float, pbar: float,
         E_U0: numpy.ndarray, V_U0: numpy.ndarray, **kwargs) → None
```

Parameters

- **R** (*float*) – ratio of amplification probabilities, R
- **pbar** (*float*) – geometric mean of amplification probabilities, \bar{p}
- **E_U0** (*np.ndarray*) – initial expected values of both strands $\mathbb{E}[\mathbf{U}_0]$, 2 by 1 matrix
- **V_U0** (*np.ndarray*) – initial variances of both strands $\text{Var}[\mathbf{U}_0]$, 2 by 2 matrix
- **f_plus** (*np.ndarray*) – molar fluorescences for each cycle/well of active reporter \mathbf{f}^+ , n by number of wells matrix Determined in units of fluorescence divided by (pmol/L)
- **f_minus** (*np.ndarray*) – molar fluorescences for each cycle/well of inactive reporter \mathbf{f}^- , n by number of wells matrix Determined in units of fluorescence divided by (pmol/L)
- **C** (*float*) – Total concentration of reporters in pmol/L, C
- **Vol** (*float*) – Total volume of solution in L, V
- **kwargs** (*dict*) – extra kwargs for Amplification class

calculate()

Performs calculations, filling in E_F, V_F, E_DX, and V_DX.

get_Cov_XiX0 (*i*)

Parameters *i* (*int*) – cycle number

Returns $\text{Cov}[X_i, X_0]$

Return type float

get_E_DX (*i*) → float

Parameters *i* (*int*) – cycle number

Returns $\mathbb{E}[\Delta X_i]$

Return type float

get_V_DX (*i*)

Parameters *i* (*int*) – cycle number

Returns $\text{Var}[\Delta X_i]$

Return type float

`src.kinetic_PCR.plot_fluorescence_curve` (*kls*: `src.kinetic_PCR.HydrolysisProbes`, *ax*,
well='A1')

Parameters

- **kls** (`HydrolysisProbes`) – Contains all amplification data. Has already computed values
- **ax** (`matplotlib.axes`) – axes to plot on
- **well** (`str`) – name of well to plot

Notes

- Plots expected value of fluorescence and shades above and below with 3 standard deviations
- Also plots maximum expected value and minimum expected value of all wells in blue dashed lines

PLOTting FIGURES

```
plot_figures.plot_figure2 (R=0.9, pbar=0.85, max_cycle=4)
```

Parameters

- **R** (*float, optional*) – choice for value of R , defaults to 0.9
- **pbar** (*float, optional*) – choice for value of \bar{p} , defaults to 0.85
- **max_cycle** (*int, optional*) – choice for maximum cycle to plot, defaults to 4

```
plot_figures.plot_figure4 (plus: src.molar_fluorescence.MolarFluorescence, minus:
                             src.molar_fluorescence.MolarFluorescence)
```

Plot figure 4

Parameters

- **plus** (*MolarFluorescence*) – fluorescence data associated with active probe
- **minus** (*MolarFluorescence*) – fluorescence data associated with inactive probe

```
plot_figures.plot_figure5 (f_plus: numpy.ndarray, f_minus: numpy.ndarray, C=0.125, Vol=2e-05,
                           R=1.0, pbar=0.9)
```

Parameters

- **f_plus** (*np.ndarray*) – molar fluorescences for each cycle/well of active reporter f^+ , n by number of wells matrix Determined in units of fluorescence divided by (pmol/L)
- **f_minus** (*np.ndarray*) – molar fluorescences for each cycle/well of inactive reporter f^- , n by number of wells matrix Determined in units of fluorescence divided by (pmol/L)
- **C** (*float, optional*) – concentration in pmol/L, defaults to 0.125
- **Vol** (*float, optional*) – volume in L, \mathcal{V} , defaults to 20e-6
- **R** (*float, optional*) – square root of ratio of probabilities, R , defaults to 1.0
- **pbar** (*float, optional*) – geometric mean of probabilities, \bar{p} , defaults to 0.9

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

p

`plot_figures`, [13](#)

s

`src.amplification`, [5](#)

`src.get_data`, [1](#)

`src.kinetic_PCR`, [11](#)

`src.molar_fluorescence`, [9](#)

`src.wells`, [3](#)